



The smart choice for
your wound care center

Solventum™ Snap™ Therapy System

The Solventum™ Snap™ Therapy System

Designed with wound care centers in mind

When you need disposable NPWT, the Snap Therapy System is the right choice.

- May help promote wound healing by removing excess exudate infectious material and tissue debris.
- Medicare coverage begins on Day 1, allowing earlier wound bed preparation without waiting four weeks to initiate therapy.*
- Compatible with other modalities like off-loading, skin substitutes grafts, and HBOT†, etc.



Delivers the proven benefits¹ of NPWT.

Clinical evidence

Armstrong¹

In a multicenter RCT, 132 patients with lower extremity diabetic and venous wounds were enrolled in the study. 118 patients were treated either with Snap Therapy System (n=59) or Solventum™ V.A.C.® Therapy (n=56), with 115 patients completing the study.

- Patients were treated for up to 16 weeks or complete wound closure.
- Primary end point analysis of wound size reduction found Snap Therapy System treated subjects demonstrated non-inferiority to V.A.C.® Therapy subjects at 4, 8, 12 and 16 weeks (p=0.0030, 0.0130, 0.0051 and 0.0044, respectively).
- The study indicated that the effect of the Snap Therapy System was not significantly different than that of the V.A.C.® Therapy System in promoting complete wound closure in the population studied (p=0.9620).
- Snap Therapy System patients reported less interruption of activities on daily living compared to V.A.C.® Therapy patients. However, pain associated with treatment was not significantly different between treatment groups.
- Other benefits noted by the authors were shorter time to dressing application and ease of use.
- However, despite randomization, the initial wound size was significantly greater in the V.A.C.® Therapy patients than in the Snap Therapy System patients (mean of 9.95cm² vs 5.37cm²; p=0.0093).

Marston²

In a multicenter RCT, 40 patients with venous leg ulcers were treated either with Snap Therapy System (n=19) or V.A.C.® Therapy (n=21).

- Patients were evaluated for 16 weeks or complete wound closure.
- Primary end point analysis of wound size reduction found Snap Therapy System treated subjects significantly greater wound size reduction than in V.A.C.® Therapy subjects at 4, 8, 12 and 16 weeks (p=0.0039, 0.0086, 0.0002, and 0.0005, respectively).
- 53% of Snap Therapy System patients achieved 50% wound closure at 30 days compared to 24% of V.A.C.® Therapy patients.
- However, despite randomization, the initial wound size was significantly greater in the V.A.C.® Therapy patients than in the Snap Therapy System patients (mean of 11.6cm² vs 4.49cm²).

*May vary by payer

†Dressings only



Easy to use for the clinician

Designed to simplify advanced wound therapy.

- Application times <10 minutes¹
- When used with the Solventum™ Snap™ Bridge Dressing Kit with Hydrocolloid Skin Barrier, provides fast and easy sealing on uneven skin surfaces and challenging body contours.
- Off-the-shelf availability



Reduced overall cost of care versus modern dressings³

The only mechanically powered NPWT system reimbursed in Wound Care Clinics.**

- Single use, disposable system

** CPT® Code 97607/8



Preserves patient quality of life¹

Negative pressure therapy that won't slow patients down.

- Minimal interference with overall activity, sleep and social interactions.¹
- Unrestricted mobility enables patients to move freely while managing their wounds.
- Small, silent, discreet design can be hidden under clothing.

Venous leg ulcer



Start of Solventum™ Snap™ Therapy






Week 4








Week 7

Diabetic hallux wound



Ordering information

	SKU	Description	Size	Quantity
	SNPA125US	Solventum™ Snap™ Therapy Cartridge -125 mmHg	60 mL	Eaches
	SNPA125US/10	Solventum™ Snap™ Therapy Cartridge -125 mmHg	60 mL	Case of 10
	BKTF14X11	Solventum™ Snap™ Bridge Dressing Kit, Foam	14 cm x 11 cm	Eaches
	BKTF14X11/10	Solventum™ Snap™ Bridge Dressing Kit, Foam	14 cm x 11 cm	Case of 10
	BKTF14X11S	Solventum™ Snap™ Bridge Dressing Kit with Hydrocolloid Skin Barrier	14 cm x 11 cm	Eaches
	BKTF14X11S/10	Solventum™ Snap™ Bridge Dressing Kit with Hydrocolloid Skin Barrier	14 cm x 11 cm	Case of 10
	SKTF10X10	Solventum™ Snap™ Advanced Dressing Kit, Foam	10 cm x 10 cm	Eaches
	SKTF10X10/10	Solventum™ Snap™ Advanced Dressing Kit, Foam	10 cm x 10 cm	Case of 10
	SKTF15X15	Solventum™ Snap™ Advanced Dressing Kit, Foam	15 cm x 15 cm	Eaches
	SKTF15X15/10	Solventum™ Snap™ Advanced Dressing Kit, Foam	15 cm x 15 cm	Case of 10
	STPAS	Solventum™ Snap™ Therapy Strap, Small	18"	Eaches
	STPAM	Solventum™ Snap™ Therapy Strap, Medium	21"	Eaches
	STPAL	Solventum™ Snap™ Therapy Strap, Large	24"	Eaches
	SRNG10	Solventum™ Snap™ Hydrocolloid Skin Barrier	2" diameter	Case of 10

For additional information, please contact your local representative.

References

1. Armstrong DG, Marston WA, Reyzelman AM, et al. Comparative effectiveness of mechanically and electrically powered negative pressure wound therapy devices: a multicenter randomized controlled trial. *Wound Rep Reg* 2012; 20(3):332-341.
2. Marston WA, Armstrong DG, Reyzelman AM, Kirsner RS. A Multicenter Randomized Controlled Trial Comparing Treatment of Venous Leg Ulcers Using Mechanically Versus Electrically Powered Negative Pressure Wound Therapy. *Advances in Wound Care*. 2015;4(2):75-82. doi:10.1089/wound.2014.0575.
3. Hutton DW, Sheehan P. Comparative effectiveness of the Snap wound care system. *Int Wound J*. 2011 Apr;8:196-205.

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.



Solventum Advanced Wound Care
12930 IH10 W
San Antonio, TX 78249
USA

Phone 800-275-4524
Web Solventum.com

© Solventum 2025. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates.

70-2013-1021-9