




Start smart with Solventum™ Veraflo™ Therapy

Backed by meta-analysis clinical
and health economic data



Key reasons to start with Veraflo Therapy: The sooner, the better¹

A meta-analysis of 13 comparative studies comprising 720 patients with a wide variety of wound types showed a positive effect when using NPWTi-d versus the control.¹



Global panel recommendations

Based on the 2020 NPWT with instillation consensus document by Paul Kim, DPM, MS et al. in the publication International Wound Journal, Veraflo Therapy should be considered early in treatment based on data and experience.³

Appropriate wound characteristics:

- Clean wounds
- Adequately cleansed and debrided wounds
- Wounds with heavy bioburden
- Contaminated/infected wounds
- Wounds that are difficult to granulate



Recommended therapy settings:

- NPWT Time: 2-3 hours

Dwell time:	NPWT setting:	Solution:
10 minutes	-125 mmHg	Saline (or other compatible solutions)


- Dressing: Solventum™ Veraflo™ Dressing unless non-viable tissue is present. Consider Solventum™ Veraflo Cleanse Choice™ Dressing as a wound cleansing option before surgical debridement or if debridement is delayed or not possible

Discontinue therapy when clinical goals are met.

For potentially significant savings



Based upon the meta-analysis by Allen Gabriel, MD et al. an economic model¹ was developed to compare the cost using Veraflo Therapy to traditional wound care options including Solventum™ V.A.C.® Therapy.



Improved outcomes can result in health economic benefits



Solventum innovation advancing NPWT because there is more than NPWT for your patients

Recent technology upgrade to the Solventum™ V.A.C.® Ulta™ Therapy Unit provides:

- Solventum™ Smart Instill™ Feature automates steps and delivers easier and less time-consuming initiation of therapy
- Automatic troubleshooting guidance
- Instill phase postpone feature
- Therapy inactive alarm timer to provide dressing changes without alarms

*These results were shown in a porcine study. It is uncertain how porcine data may correlate to human results.

Solventum™ V.A.C. Ulta™ Therapy System ordering information for Veraflo Therapy

Pump & components		Kits & dressings	
ULTDEV01/US	V.A.C.® Ulta Therapy Unit, United States	ULTVFL05SM	Solventum™ Veraflo™ Dressing, small, 5-pack
ULTLNK0500	Solventum™ Instillation Cassette, 5-pack	ULTVFL05MD	Solventum™ Veraflo™ Dressing, medium, 5-pack
ULTDUO0500	Solventum™ VeraT.R.A.C. Duo™ Tube Set, 5-pack	ULTVFL05LG	Solventum™ Veraflo™ Dressing, large, 5-pack
M8275063/5	V.A.C.® Ulta Therapy Unit – Solventum™ V.A.C.® Canister with Gel, 500 mL	DTAC10LDP	Solventum™ Dermatac™ Drape, Case of 10
M8275093/5	V.A.C.® Ulta Therapy Unit – Solventum™ V.A.C.® Canister with Gel, 1000 mL		

Start smart with Veraflo Therapy, visit [Solventum.com/Veraflo](https://www.solventum.com/Veraflo)

References:

1. Gabriel A, Camardo M, O'Rourke E, Gold R, Kim P J. Effects of Negative-Pressure Wound Therapy With Instillation versus Standard of Care in Multiple Wound Types: Systematic Literature Review and Meta-Analysis. *Plast Reconstr Surg.* 2021 Jan 1;147(1S-1):68S-76S. doi: 10.1097/PRS.0000000000007614. PMID: 33347065.
2. Kim PJ, Lookess S, Bongards C, Griffin LP, Gabriel A. Economic model to estimate cost of negative pressure wound therapy with instillation vs control therapies for hospitalised patients in the United States, Germany, and United Kingdom. *Int Wound J.* 2021;1-7. doi:10.1111/iwj.13689.
3. Kim PJ, Attinger CE, Constantine T, et al. Negative pressure wound therapy with instillation: international consensus guidelines update. *Int Wound J.* 2020 Feb;17(1):174-186.
4. Carroll C A, Kilpadi D, Ingram S. Comparison of a Dressing with Through Holes to Alternative Reticulated Open Cell Foam Dressings in a Porcine Model with Slough and Treatment with Negative Pressure Wound Therapy with Instillation. Poster presented at the Symposium on Advanced Wound Care, October 20-22, 2017; Las Vegas, NV.
5. Solventum™ V.A.C.® Ulta Negative Pressure Wound Therapy System Safety Information and Solventum™ V.A.C. Veraflo Cleanse Choice™ Dressing System Application Instructions.



Solventum Advanced Wound Care

12930 IH10 W
San Antonio, TX 78249

Phone 800-275-4524

Web [Solventum.com](https://www.solventum.com)

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals. Rx only.

© Solventum 2026. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates.

70-2011-8181-8