

3M™ Tube Securement Device

This strong, easy-to-use solution secures non-critical medical tubes made from silicone, PVC and silicone coated latex spanning 6–24 French. The tube securement device provides additional securement to fluid drainage tubing to help stabilise and prevent movement. It is gentle to skin and designed to support pressure injury and medical adhesive-related skin injury (MARSI) prevention practices.

Sterile device

In single-use packaging.

Purposely shaped


Designed to lift tubes off skin and help promote consistent application.

Reinforced soft cloth backing

Provides extra strength and flexibility.

Two sizes

Small device for medical tubes and catheters 6–12 Fr and large device for medical tubes and catheters 12–24 Fr.



Skin performance: The feather designates products that deliver the securement power you need while minimising damage to skin.

Reliable securement

Unique adhesive secures medical tubes and catheters for up to 4 days.¹

Ease of use

Easy to apply and removes cleanly in one piece; allows for standardisation of technique.

Skin safe

Non-sensitising silicone adhesive has low potential to cause dermatitis with minimal skin damage and pain upon removal.¹

Designed to help reduce pressure injury risk

With no hard components to cause pressure, product's unique shape minimises contact between skin and tube.

¹ 3M Data on File, EM-05-014960.

Considerations and indications

When it comes to securing tubing, there are a variety of factors to keep in mind, including securement performance, skin protection and clinician ease of use. The 3M™ Tube Securement Device is indicated for non-critical tube securement or in conjunction with sutures, primary securement dressings and/or urinary catheter balloons per facility protocol. Two sizes accommodate tubes spanning 6-24 French.

Applications



Vascular access

The tube securement device meets the Infusion Therapy Standards of Practice guidelines as an adhesive securement device and should be used in conjunction with sutures and/or primary securement dressings.

Consider using the securement device as secondary securement for primary tubing and extension tubing on central lines and peripheral lines. It can help organise lines with catheters that have multiple lumens.



Surgical drains

Provides additional securement in conjunction with sutures to fluid drainage tubing to help prevent movement.



Urinary catheters



When compared to StatLock® Foley Stabilisation Device in a comparative prospective study of 24 healthy human subjects showed that the tube securement device:

- 9x faster removal. The average removal time for the tube securement device was 13 seconds compared to 126 seconds for StatLock® Foley Stabilisation Device (P<.0001).*
- No removal agent. On average no wipes were needed to remove the tube securement device, while the StatLock® Foley Stabilisation Device requires 12 alcohol wipes for removal.*

Feel secure in your choice. Trust the science behind tube securement device for reliable performance.

Ordering information

Refer to product Instructions for Use (IFU) for other important information.

Product	Product name	Product number	Tube & catheter size	Product size	Devices/box	Boxes/case
	3M™ Tube Securement Device, Small	1500U	6-12 French	4.3 cm x 6.1 cm	25	4
	3M™ Tube Securement Device, Large	1501U	12-24 French	5.5 cm x 9.0 cm	25	4

*As tested in a comparative, prospective clinical study on 24 healthy human subjects after 3 days of wear followed by device removal from skin per IFU.



KCI Medical Australia Pty Ltd
Level 3, Building A, 1 Rivett Road, North Ryde 2113
Phone: 1800 945 183
KCI New Zealand Unlimited
Suite 1701, Level 17 PwC Tower, 15 Customs St West Auckland 1010
Phone: 09 884 2001

To order, contact your
Solventum Representative.

© Solventum 2025. Solventum and the S logo are trademarks of Solventum or its affiliates. 3M and the 3M logo are licensed trademarks of 3M. All rights reserved.