



Manage small, low-
exudating wounds with
disposable NPWT¹

Solventum[™] Snap[™] Therapy System





Solventum™ Snap™ Therapy System

Combines the simplicity of advanced wound dressings with the proven benefits¹ of negative pressure wound therapy in a discreet design.

Clinical evidence

Armstrong¹

In a multicenter RCT, 132 patients with lower extremity diabetic and venous wounds were enrolled in the study. 118 patients were treated either with Snap therapy system (n=59) or Solventum™ V.A.C.® Therapy Unit (n=56), with 115 patients completing the study.

- Patients were treated for up to 16 weeks or complete wound closure.
- Primary end point analysis of wound size reduction found that Snap therapy system treated subjects demonstrated non-inferiority to V.A.C.® Therapy subjects at 4, 8, 12 and 16 weeks ($p=0.0030$, 0.0130 , 0.0051 and 0.0044 , respectively).
- The study indicated that the effect of the Snap therapy system was not significantly different than that of the V.A.C. Therapy System in promoting complete wound closure in the population studied ($p=0.9620$).
- Snap therapy system patients reported less interruption of activities on daily living compared to V.A.C. Therapy patients. However, pain associated with treatment was not significantly different between treatment groups.
- Other benefits noted by the authors were shorter time to dressing application and ease of use.

- However, despite randomization, the initial wound size was significantly greater in the V.A.C. Therapy patients than in the Snap therapy system (mean of 9.95cm^2 vs 5.37cm^2 ; $p=0.0093$)

Marston³

In a multicenter RCT, 40 patients with venous leg ulcers were treated either with Snap therapy system (n=19) or V.A.C. Therapy (n=21)

- Patients were evaluated for 16 weeks or complete wound closure.
- Primary end point analysis of wound size reduction found that Snap therapy system treated subjects had significantly greater wound size reduction than in V.A.C. Therapy subjects at 4, 8, 12 and 16 weeks ($p\text{-value}=0.0039$, 0.0086 , 0.0002 , and 0.0005 , respectively).
- 53% of Snap therapy system patients achieved 50% wound closure at 30 days compared to 24% of V.A.C. Therapy patients.
- However, despite randomization, the initial wound size was significantly greater in the V.A.C. Therapy patients than in the Snap therapy system patients (mean of 11.6cm^2 vs 4.49cm^2).

Solventum™ Snap™ Therapy System

Silent cartridge technology

All the power of negative pressure wound therapy without the help of a battery.



Patient-friendly features


- Small, lightweight design can be hidden under clothes
- No batteries to interfere with daily living
- Visual indicator displays when cartridge is full or discharged

Innovative design

- Proprietary spring mechanism generates consistent, even levels of pressure²
- Technology gels exudate for exudate containment
- -125 mmHg pressure setting

Solventum™ Snap™ Advanced Dressing Kits

Proprietary hydrocolloid dressing offers periwound protection and easy removal. Additional accessories are designed to simplify dressing applications.




Snap Advanced Dressing Kits

- Absorptive hydrocolloid properties help to maintain seal in the presence of exudate or sweat to help reduce periwound maceration
- Disposable components with off-the-shelf convenience
- Fully-integrated microport enables flexibility and a tight bending radius for wounds located in difficult areas
- Cut-to-length tubing and integrated one-way flow for improved safety
- Available in two sizes:
10 cm x 10 cm and 15 cm x 15 cm

Solventum™ Snap™ Bridge Dressing Kit

- Completely flat, comfortable dressing surface to help minimize further pressure damage
- Integrated bridge and port for one-step application
- Soft pad cushioning under bridge for improved patient comfort

Usage Locations:



Interface layers

- Available with a blue foam interface
- Facilitates even levels of negative pressure²



Solventum™ Snap™ Bridge Dressing Kit with Hydrocolloid Skin Barrier

- Fast and easy sealing on uneven skin surfaces and challenging body contours
- Reduces accessories needed to seal and protect the wound from moisture
- Increased adhesion of the Solventum™ Snap™ Dressing on dry and uneven skin



Solventum™ Snap™ Therapy Strap

- Soft strap enables device to be conveniently and discreetly worn under clothes

Solventum™ Snap™ Therapy System ordering information

Solventum™ Snap™ Therapy Cartridge			
Catalog number	Pressure	Capacity	Quantity
SNPA125	-125 mmHg	60 mL	1 each
SNPA125/10	-125 mmHg	60 mL	10 each

Solventum™ Snap™ Advanced Dressing Kit			
Catalog number	Size	Interface	Quantity
SKTF10X10	10 cm x 10 cm	Foam	1 each
SKTF10X10/10	10 cm x 10 cm	Foam	10 each
SKTF15X15	15 cm x 15 cm	Foam	1 each
SKTF15X15/10	15 cm x 15 cm	Foam	10 each

Solventum™ Snap™ Bridge Dressing Kit			
Catalog number	Size	Interface	Quantity
BKTF14X11	14 cm x 11 cm	Foam	1 each
BKTF14X11/10	14 cm x 11 cm	Foam	10 each
BKTF14X11S	14 cm x 11 cm with Solventum™ Snap™ Bridge Dressing Kit with Hydrocolloid Skin Barrier	Foam	1 each
BKTF14X11S/10	14 cm x 11 cm with Snap Bridge Dressing Kit with Hydrocolloid Skin Barrier	Foam	10 each

Solventum™ Snap™ Hydrocolloid Skin Barrier		
Catalog number	Size	Quantity
SRNG10	2" diameter	10 each

Solventum™ Snap™ Therapy Strap		
Catalog number	Size	Quantity
STPAS	Small 18"	1 each
STPAM	Medium 21"	1 each
STPAL	Large 24"	1 each

To order product or for more information, please contact your local representative.

References

1. Armstrong DG, Marston WA, Reyzelman AM, Kirsner RS. Comparative effectiveness of mechanically and electrically powered negative pressure wound therapy devices: a multicenter randomized controlled trial. *Wound Rep Reg* 2012; 20(3):332-341. 2. Fong KD, Hu D, Eichstadt S et al. The SNaP system: biomechanical and animal model testing of a novel ultraportable negative-pressure wound therapy system. *Plastic and Reconstructive Surgery*. 2010 May;125(5):1362-71. 3. Marston WA, Armstrong DG, Reyzelman AM, Kirsner RS. A Multicenter Randomized Controlled Trial Comparing Treatment of Venous Leg Ulcers Using Mechanically Versus Electrically Powered Negative Pressure Wound Therapy. *Advances in Wound Care*. 2015;4(2):75-82. doi:10.1089/wound.2014.0575.

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.



Solventum Advanced Wound Care
12930 IH10 W
San Antonio, TX 78249
USA

Phone 1-800-275-4524
Web Solventum.com

© Solventum 2026. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates. Other trademarks are the property of their respective owners.

70-2011-8245-1 PRA-PM-US-00289 (12/21)