


































Product	Code	Packaging	Indications for use	Time to result	Incubator/auto-reader
Super rapid readout	 <b>1491</b> Dark blue cap	50 each/box 4 boxes/case	 <b>Gravity:</b> 132°C = 3 or 10 minutes 135°C = 3 or 10 minutes		 490/490M Auto-readers
	 <b>1492V</b> Brown cap	50 each/box 4 boxes/case	 <b>Dynamic-air-removal (pre-vacuum and SFPP):</b> 132°C = 3, 3.5, 4, 5.5 or 6 minutes 134°C = 3 or 4 minutes 135°C = 3, 3.5 or 10 minutes		
	 <b>1493</b> Light blue cap	50 each/box 4 boxes/case	 <b>Gravity:</b> 121°C = 30 minutes <b>Dynamic-air-removal (pre-vacuum and SFPP):</b> 121°C = 15, 20, 30 or 35 minutes		
	 <b>1493P</b> Light blue cap	20 each/box 5 boxes/case			
	 <b>1492PCD-24</b>	24 packs 24 controls/case	 <b>Dynamic-air-removal (pre-vacuum and SFPP):</b> 132°C = 4 minutes 134°C = 4 minutes 135°C = 3 minutes		
	 <b>1492PCD-05</b>	24 packs 5 controls/case			
 <b>1493PCDG-24</b>	24 packs 24 controls/case	 <b>Gravity displacement:</b> 121°C = 30 minutes			
Rapid readout	 <b>1291</b>	50 each/box 4 boxes/case	 <b>Gravity:</b> 132°C ≥ 3 minutes		 390 Auto-reader
	 <b>1292</b>	50 each/box 4 boxes/case	 <b>Gravity:</b> 121°C ≥ 20 minutes <b>Vacuum-assisted:</b> 132°C ≥ 4 minutes		
	 <b>1296</b>	25 packs/box 25 controls/box	 <b>Gravity:</b> 121°C ≥ 40 minutes <b>Vacuum-assisted:</b> 132°C ≥ 4 minutes		
Standard	 <b>1261</b>	100 each/box 4 boxes/case	 <b>Gravity:</b> 132°C		 118 Incubator
	 <b>1262</b>	100 each/box 4 boxes/case	 <b>Gravity:</b> 121°C <b>Vacuum-assisted:</b> 132°C <b>Gravity (table-top sterilisers only):</b> 132°C > 3 minutes		
	 <b>1262P</b>	25 each/box 4 boxes/case			

\*3M™ Attest™ Biological Indicators 1491, 1492V, 1493 can be read in 24 minutes when used with 3M™ Attest™ Mini Auto-reader 490M or with 3M™ Attest™ Auto-reader 490 or 490H with software version 4.0.0 or greater.

\*\*Must be used in conjunction with 3M™ Attest™ Biological Indicator Incubator 118.