

	Order code	Dispensing system	Viscosity low  high	Impression technique	Recommended max. working time (min:sec)	Intra-oral setting time at 37°C (min:sec)
Super quick setting materials						
3M™ Impregum™ Penta™ Super Quick (Heavy Body)	69406			1-step	0:45	2:00
3M™ Impregum™ Penta™ Super Quick (Medium Body)	69385			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Super Quick (Medium Body)	69417			Monophase or 1-step	0:45	2:00
3M™ Impregum™ Super Quick (Light Body)	69379			1-step	0:45	2:00
3M™ Impregum™ Super Quick (Ultra Light Body)	69418			1-step	0:45	2:00
Quick setting materials						
3M™ Impregum™ Penta™ Soft Quick (Heavy Body)*	31769			1-step	1:00	3:00
3M™ Impregum™ Penta™ Soft Quick (Medium Body)	31765			Monophase or 1-step	1:00	3:00
3M™ Impregum™ Garant™ Soft Quick (Medium Body)	31788			Monophase or 1-step	1:00	3:00
3M™ Impregum™ Soft Quick (Light Body)	31767			1-step	1:00	3:00
Regular setting materials						
3M™ Impregum™ Penta™ H DuoSoft™ (Heavy Body)	31740			1-step	2:30	3:30
3M™ Impregum™ Penta™ Soft (Medium Body)	31730			Monophase or 1-step	2:45	3:15
3M™ Impregum™ Garant™ Soft (Medium Body)	31785			Monophase or 1-step	1:45	4:15
3M™ Impregum™ Penta™ (Medium Body)	31644			Monophase	2:45	3:15
3M™ Impregum™ Penta™ L DuoSoft™ (Light Body)	31745			1-step	3:15	3:15
3M™ Impregum™ Garant™ L DuoSoft™ (Light Body)	31751			1-step	2:00	3:30

*3M™ Impregum™ Penta™ Soft Quick Heavy Body (31769) is not available in New Zealand.

Always follow the Instructions for Use (IFU) and refer to IFU for full indications, precautions and warnings.

© Solventum 2026. Solventum, the S logo, Garant, Impregum, Imprint and Penta are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. All other trademarks are owned by their respective owners.

	Order code	Dispensing system	Viscosity low  high	Max. working time (23°C) 	Max. intra-oral syringe time (37°C) 	Intra-oral setting time (37°C) 
Super quick setting materials						
3M™ Imprint™ 4 Penta™ Super Quick Heavy	71485			1:15	–	1:15
3M™ Imprint™ 4 Super Quick Heavy	71493			1:15	–	1:15
3M™ Imprint™ 4 Super Quick Ultra-Light	71491			–	0:35	1:15
3M™ Imprint™ 4 Super Quick Light	71490			–	0:35	1:15
3M™ Imprint™ 4 Super Quick Regular	71489			–	0:35	1:15
Regular setting materials						
3M™ Imprint™ 4 Penta™ Putty	71486			1:30	–	2:30
3M™ Imprint™ 4 Penta™ Heavy	71484			2:00	–	2:00
3M™ Imprint™ 4 Heavy	71492			2:00	–	2:00
3M™ Imprint™ 4 Light	71488			–	1:00	2:00
3M™ Imprint™ 4 Regular	71487			–	1:00	2:00
Alginate replacement						
3M™ Imprint™ 4 Penta™ Preliminary Super Quick	71522			1:00	–	1:30
3M™ Imprint™ Preliminary Super Quick	71540			1:00	–	1:30
3M™ Imprint™ 4 Penta™ Preliminary	71521			1:30	–	2:00
3M™ Imprint™ Preliminary*	71539			1:30	–	2:00

*3M™ Imprint™ Preliminary (71539) is not available in New Zealand.