


3M™ RelyX™ Universal Resin Cement &
3M™ Scotchbond™ Universal Plus Adhesive

The two that can do it all

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement – and as an adhesive cement together with 3M™ Scotchbond™ Universal Plus Adhesive. The adhesive serves as a universal primer for all materials and covers all direct and indirect bonding indications.

Eliminate the hassle and confusion of multiple resin cements, primers and adhesives and simplify your direct and indirect restorative workflows. Less products to keep on stock and replenish. Saves space, time and money.

 First radiopaque all-in one universal adhesive



 Further boosts bond strength of RelyX Universal Resin Cement

 Fully aligned system: adhesive cured by cement, no light-cure needed

 Universal primer for all restorative materials



Game-changing syringe Proprietary 3M design

- Micro mixing tip providing 80% less cement waste*

- 15 - 20 applications from 3.4 g**

- Hygienic self-sealing syringe – for storage without the used mixing tip



- More ergonomic shape with 50% less plastic waste*

- Thin, long elongation tip for easy root canal application




- Easy excess clean-up



Image courtesy of Dr. Jan Gueth, Germany

- Excellent self-adhesive bond strength to dentine in self- and light-cure

Shear bond strength (MPa) to dentin after artificial aging



R. Afutu, M. Abreu, G. Kugel, Tufts University, *J. Dent. Res.* Vol 98A, No 3629, 2019

* Per application compared to currently available standard automix systems.

Source: M. Cowen, J.M. Powers, Dental Advisor, Biomaterials Research Report, Number of Automix Applications and Mixing Efficiency, No. 134 – June 23, 2020.

** Based upon the amount of cement used in a restoration from 60 µL (e.g. inlay, onlay, veneer) to 90 µL (e.g. crown).

High aesthetics with four colour-stable, fluorescent shades



Ordering information

3M™ RelyX™ Universal Resin Cement

Item no.	Product information		
Intro kit			
56968	4 x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4 g each), 60 x 3M™ RelyX™ Universal Micro Mixing Tips, 20 x elongation tips, 1 x vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 x orange microbrushes, step-by-step guides		
Trial kits			
Each trial kit contains: 1 x 3M™ RelyX™ Universal Resin Cement syringe (3.4 g), 15 x 3M™ RelyX™ Universal Micro Mixing Tips, 5 x elongation tips, 1 x vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 x blue etchant tips, 50 x orange microbrushes, step-by-step guides			
56969	Shade TR	56970	Shade A1
Refills			
Each refill contains: 1 x 3M™ RelyX™ Universal Resin Cement syringe (3.4 g), 15 x 3M™ RelyX™ Universal Micro Mixing Tips, 5 x elongation tips			
56971	Shade TR	56973	Shade A3O
56972	Shade A1	56974	Shade WO
Value packs			
Each value pack contains: 3 x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45 x 3M™ RelyX™ Universal Micro Mixing Tips, 15 x elongation tips			
56977	Shade TR	56978	Shade A1
Mixing tips			
56975	30 x 3M™ RelyX™ Universal Micro Mixing Tips		
56976	15 x 3M™ RelyX™ Universal Micro Mixing Tips with elongation tips		



3M™ RelyX™ Universal Resin Cement Intro Kit (56968)



3M™ RelyX™ Universal Resin Cement Trial Kit TR (56969)



3M™ RelyX™ Universal Resin Cement Refill TR (56971)



3M™ Scotchbond™ Universal Plus Adhesive Refill (41294)

3M™ Scotchbond™ Universal Plus Adhesive

Vial			
41293	Intro kit: 1 x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 x disposable brushes, 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 x dispensing tips		
Refill: 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)			
41294	1 x 5 ml vial	41295	3 x 5 ml vial
Unit dose			
3M™ Scotchbond™ Universal Plus Adhesive Unit Dose L-Pops, 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 x dispensing tips			
41297	Intro kit: 50 x L-Pop	41299	200 x L-Pop
41298	100 x L-Pop		

3M™ Scotchbond™ Universal Etchant Etching Gel

41263	2 x 3M™ Scotchbond™ Universal Etchant (3 ml each), 50 x dispensing tips		
-------	---	--	--



Solvantum Dental Solutions

KCI Medical Australia Pty Ltd
Level 3, Building A, 1 Rivett North Ryde, NSW 2113
1300 363 454

KCI New Zealand Unlimited
Suite 1701, Level 17, PwC Tower, 15 Customs Street West
Auckland Central, Auckland 1010
0800 808 182

www.solvantum.com

Always follow the Instructions for Use (IFU) and refer to IFU for full indications, precautions and warnings. Claims supported by 3M data on file.

© Solvatum 2026. Solvatum, the S logo and RelyX are trademarks of Solvatum or its affiliates. 3M, the 3M logo and Scotchbond are trademarks of 3M. All other trademarks are owned by their respective owners.