



3M™ RelyX™ Universal Resin Cement & 3M™ Scotchbond™ Universal Plus Adhesive

The two that can do it all

RelyX Universal Resin Cement works both as a standalone, self-adhesive cement – and as an adhesive cement together with Scotchbond Universal Plus Adhesive. The adhesive serves as a universal primer for all materials and covers all direct and indirect bonding indications.

Eliminate the hassle and confusion of multiple resin cements, primers and adhesives and simplify your direct and indirect restorative workflows. Fewer products to keep in stock and replenish. Saves space, time and money.

- ✓ First radiopaque all-in-one universal adhesive
- ✓ Further boosts bond strength of RelyX Universal Resin Cement
- ✓ Fully aligned system: adhesive cured by cement, no light-cure needed
- ✓ Universal primer for all restorative materials



Game-changing syringe. Proprietary Solventum design.

- ✓ Micro Mixing Tip providing 80% less cement waste*
- ✓ 15 applications from 3.4g**
- ✓ More hygienic, easier to clean***
- ✓ More ergonomic shape with 50% less plastic waste*
- ✓ Thin, long elongation tip for easy root canal application




- ✓ Easy excess clean-up



Image courtesy of Dr. Jan Gueth, Germany

- ✓ Superior self-adhesive bond strength to dentin in self- and light-cure




* per application compared to currently available standard automix systems

** on average

*** compared to currently used automix syringes

High esthetics with four color-stable, fluorescent shades



Translucent (TR) A1 A3 Opaque (A3O) White Opaque (WO)

Ordering information

3M™ RelyX™ Universal Resin Cement

Item no.	Product information		
Intro Kit			
56968	4x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60x 3M™ RelyX™ Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides		
Trial Kits			
56969	Shade TR	56970	Shade A1
Each trial kit contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides			
Refills			
56971	Shade TR	56973	Shade A3O
56972	Shade A1	56974	Shade WO
Each refill contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips			
Value Packs			
56977	Shade TR	56978	Shade A1
Each value pack contains: 3x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45x 3M™ RelyX™ Universal Micro Mixing Tips, 15 Elongation Tips			
Mixing Tips			
56975	30x 3M™ RelyX™ Universal Micro Mixing Tips		
56976	15x 3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips		



RelyX Universal Resin Cement Intro Kit (56968)



RelyX Universal Resin Cement Trial Kit TR (56969)



RelyX Universal Resin Cement Refill TR (56971)



Scotchbond Universal Plus Adhesive Refill (41294)

3M™ Scotchbond™ Universal Plus Adhesive

Item no.	Product information		
Vial			
41293	Intro Kit: 1x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips		
Refill: 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)			
41294	1x 5 ml vial	41295	3x 5 ml vial
Unit Dose			
3M™ Scotchbond™ Universal Plus Adhesive Unit Dose, 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips			
41297	Intro Kit: 50x L-Pop	41299	200x L-Pop
41298	100x L-Pop	41304	400x L-Pop

Scotchbond Universal Plus Adhesive

Unit Dose: For efficient hygiene management. Reduces the risk of cross-contamination.



Scotchbond Universal Plus Adhesive Unit Dose 100 L-Pop (41298)

3M™ Scotchbond™ Universal Etchant Etching Gel

Item no.	Product information
41263	2x 3M™ Scotchbond™ Universal Etchant (3 ml each), 50 dispensing tips



Solventum Dental Solutions

1-800-634-2249

Solventum.com

© Solventum 2026. Solventum, the S logo and RelyX are trademarks of Solventum or its affiliates. 3M, the 3M logo and Scotchbond are trademarks of 3M. Other trademarks are the property of their respective owners.