



3M™ Scotchbond™
Universal Plus Adhesive

More than just a strong bond



3M™ Scotchbond™ Universal Plus Adhesive takes universal adhesives to the next level. It keeps all the benefits of the original 3M™ Scotchbond™ Universal Adhesive, the first truly universal adhesive, while adding even greater control and predictability – because a good adhesive should provide more than just a strong bond.

From excellent...

Built on the clinically proven and world's most researched universal adhesive

All direct and indirect indications
light-curing, self-curing or dual-curing methacrylate-based composites / cements

All etching techniques
total-etch, self-etch, and selective enamel etch

Bonds to all restorative materials
without the need for separate primers

Available in vial and unit dose
for efficient hygiene management

Virtually no post-operative sensitivity
in field evaluation:

99% of dentists reported that 3M™ Scotchbond™ Universal Plus Adhesive shows virtually no post-operative sensitivity.*

*Field evaluation with 309 dentists and over 20,000 applications.

...to exceptional

Formulation carefully optimised to achieve unique features

First radiopaque universal adhesive
reducing the risk of X-ray misdiagnosis and invasive overtreatment

Bonds and seals caries-affected dentine
maximising preservation of natural tooth structure, supporting minimally invasive preparation guidelines

Enhanced Silane Functionality that offers Gold Standard adhesion to all dental substrates
and in particular, glass ceramics

Full dual- and self-cure compatibility
no need for an additional dual-cure activator bottle

BPA free formulation
free of BPA derivatives like BisGMA




Fully aligned system with 3M™ RelyX™ Universal Resin Cement
excellent bond strength for virtually all resin cement indications. 3M Scotchbond Universal Plus Adhesive is cured when in contact with RelyX Universal Resin Cement




Added radiopacity

On average, dentists see at least one difficult to interpret X-ray per day due to radiolucency under an existing restoration. This could lead to misdiagnosis and invasive overtreatment. 3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque universal adhesive with dentine-like radiopacity, helping to overcome this issue.



Non-Radiopaque Adhesive



3M™ Scotchbond™ Universal Plus Adhesive

“Radiolucent adhesive pooling underneath a restoration doubles the risk of restoration replacement.”

Dr. John O. Burgess,
University of Alabama, USA

Intentional adhesive pooling under Class I restorations – X-rays of extracted human molars, filling material 3M™ Filtek™ One Bulk Fill Restorative. Internal Data, H. Loll *et al.*, *J Dent Res* 99 (Spec Iss A): No. 757, 2020

Ordering information

Item no. Product information

3M™ Scotchbond™ Universal Plus Adhesive – Vial

41293	Intro kit vial 1 vial 5 ml, 50 disposable brushes, 1 etchant syringe 3 ml, 25 dispensing tips
41294	Refill vial – 1 x 5 ml
41295	Refill vial – 3 x 5 ml
41263	3M™ Scotchbond™ Universal Etchant Etching Gel 2 syringes x 3 ml, 50 dispensing tips

3M™ Scotchbond™ Universal Plus Adhesive – Unit Dose

41297	Intro kit unit dose 50 unit doses, 1 etchant syringe 3 ml, 25 dispensing tips
41298	Unit dose 100 100 unit doses, 1 etchant syringe 3 ml, 25 dispensing tips
41299	Unit dose 200 200 unit doses, 1 etchant syringe 3 ml, 25 dispensing tips



Vial
with smooth and even surface – easy to wipe off excess



Unit Dose
for efficient hygiene management – reduces the risk of cross contamination



Solventum Dental Solutions

KCI Medical Australia Pty Ltd
Level 3, Building A, 1 Rivett North Ryde, NSW 2113
1300 363 454

KCI New Zealand Unlimited
Suite 1701, Level 17, PwC Tower, 15 Customs Street West
Auckland Central, Auckland 1010
0800 808 182

www.solventum.com



Scan to find out more

© Solventum 2026. Solventum, the S logo, Filtek and RelyX are trademarks of Solventum or its affiliates. 3M, the 3M Logo and Scotchbond are trademarks of 3M. All other trademarks are owned by their respective owners.