

# Sof-Lex™ Contouring, Finishing and Polishing Systems

A full line of finishing and polishing solutions ... from beginning to end

The Sof-Lex system is colour-coded from dark (coarse) to light (superfine) for an easily followed step-by-step process. Easy, pop-on mounting lets you quickly change discs or spirals—so you can move efficiently through your finishing and polishing procedure.

**Sof-Lex™ Pop-on Contouring and Polishing Disc Refills - Blue series** (bag of 85 single use discs)

9.5 mm		12.7 mm	
1981C		1982C	
1981M		1982M	
	Coarse		Coarse
	Medium		Medium
1981F		1982F	
1981SF		1982SF	
	Fine		Fine
	Super fine		Super fine




**Sof-Lex™ Pop-on Extra Thin Disc Refills - Orange series** (bag of 85 single use discs)

9.5 mm		12.7 mm	
2381C		2382C	
2381M		2382M	
	Coarse		Coarse
	Medium		Medium
2381F		2382F	
2381SF		2382SF	
	Fine		Fine
	Super fine		Super fine

**Sof-Lex™ Diamond Polishing System**

Item #	Product information
 5090	Sof-Lex™ Pre-Polishing Spirals (Beige) Refill (bag of 15 multi-use spiral wheels)
 5091	Sof-Lex™ Diamond Polishing Spirals (Pink) Refill (bag of 15 multi-use spiral wheels)
5092-IM	Sof-Lex™ Diamond Polishing System - Introductory kit Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); 1 Mandrel; Instructions for Use; Technique Card

**Sof-Lex™ Finishing and Polishing Aluminium Oxide Strips**

 1954	178 mm x 4 mm Clear coarse/Medium (150 strips)
 1954N	178 mm x 2 mm Clear coarse/Medium - Narrow (100 strips)
 1956	178 mm x 4 mm Yellow fine/Super fine (100 strips)

# Sof-Lex™ Finishing and Polishing System






Sof-Lex™ Discs



Sof-Lex™ Diamond Polishing System



Sof-Lex™ Strips



## Solventum Dental Solutions

KCI Medical Australia Pty Ltd  
Level 3, Building A, 1 Rivett North Ryde, NSW 2113  
1300 363 454

KCI New Zealand Unlimited  
Suite 1701, Level 17, PwC Tower, 15 Customs Street West  
Auckland Central, Auckland 1010  
0800 808 182

[www.solventum.com](http://www.solventum.com)

Always follow the instructions for use and refer to IFU for full indications, precautions and warnings. 3M data on file supports multi-use reference.

© Solventum 2026. Solventum, the S logo and Sof-Lex are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. All other trademarks are owned by their respective owners.