

Now our most popular flowable is easier to inject


3M™ Filtek™ Supreme Flowable Restorative is available in an ergonomic syringe designed for improved comfort and easier injection.*

Reduce bubbles and waste

Virtually no bubbles or material “run-on” during dispensing gives you better control.*

No more guessing

The blue syringe barrel indicates the remaining material volume. When your syringe is empty, only the white plunger shows.




Bendable cannula

Easy-to-bend without kinking for better access.*

Easy to hold and inject

Triangle finger plate and plunger make it quick and easy to personalize your grip.*



Key benefits

- Virtually no bubbles or material “run-on”*
- Easy-to-bend cannula for improved access*
- Ergonomic design is easy to hold and inject*
- Excellent adaptation, polish retention and wear resistance*
- Versatile for a wide range of indications and techniques

Can be warmed for injection moulding

Warm up to 70°C for up to one hour, up to 25 times.*



Clinical cases

3M™ Filtek™ Supreme Flowable Restorative offers a smooth transition from material to tooth for aesthetic restorations.

Photos courtesy of Dr. Gunnar Reich, Munich, Germany



Before

Flowable placement

After

3M™ Filtek™ Supreme Flowable Restorative adapts easily and seals the margin.

Photos courtesy of Dr. Giuseppe Chiodera, Brescia, Italy




Before

Flowable placement

After

Indications

(Incremental up to 2 mm)



Refer to Instructions for Use for full list of indications.

Ordering information



3M™ Filtek™ Supreme Flowable Restorative Refills

Syringe refills include: 2 – 2g syringes; 20 dispensing tips.

| Shade | Item no. |
|-------------|----------|
| A1 | 6032A1 |
| A2 | 6032A2 |
| A3 | 6032A3 |
| A3.5 | 6032A3.5 |
| A4 | 6032A4 |
| B1 | 6032B1 |
| B2 | 6032B2 |
| C2 | 6032C2 |
| D2 | 6032D2 |
| White | 6032W |
| Extra White | 6032XW |
| OA3 | 6032AO3 |

| Item no. | Product info. |
|----------|---|
| 3700T | 3M™ Filtek™ Supreme Flowable Restorative Dispensing Tips, 20-gauge, Gray (20)  |

References

1. Including small, non stress-bearing occlusal restorations.
2. Repair of acrylic temporary and composite resin.

*3M Internal data



Solventum
Charnwood Campus, Bakewell Road Loughborough, LE11 5TH, United Kingdom

Solventum 3M Healthcare Ireland Ltd
2 Cumberland Place, Fenian Street, Dublin 2, D02 H0V5, Ireland

For more information, visit [Solventum.com](https://www.solventum.com)

© Solventum 2025. Solventum and the S logo are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. All other trademarks are owned by their respective owners.