

3M Science.
Applied to Life.™

3M™ Clarity™ Ceramic Brackets

**Choose brackets
that give you choice.**



For bracket treatment, the choice is clear.

Available ligated or self-ligating for maximum choice, 3M™ Clarity™ Ceramic Brackets give your practice all the benefits of esthetic treatment without compromising performance.

3M™ Clarity™
Advanced Ceramic
Brackets



3M™ Clarity™
Ultra Self-Ligating
Brackets



Predictable squeeze debonding.

Reduce the uncertainty of the debonding process with a unique, proprietary stress concentrator and a simple squeeze.



Trusted durability.

Made of fine-grained alumina material, designed for strength.



Designed for esthetic comfort.

Smooth surfaces and rounded corners reduce rubbing and irritation, without sacrificing esthetics.



Easy placement.

With washable horizontal and vertical reference markers, bracket placement is more efficient than ever.

Learn more at [3M.com/CeramicBrackets](https://www.3M.com/CeramicBrackets)

3M™ APC™ Flash-Free Adhesive

Simple, efficient bonding process.

Clarity Ceramic Brackets are available precoated with APC Flash-Free Adhesive, the most efficient bonding system in orthodontics, to provide a better bonding experience to you and your patients.



Average bracket bond

failure rate of less than 2%^{**1}



40% reduction in time

required for bracket placement, positioning, and clean-up.¹



Up to 24% faster clean up

than traditional adhesive removal.¹



Discover the difference at [3M.com/APCFF](https://www.3M.com/APCFF)

^{**}Within the first ninety days of bonding. 1. 3M data on file.

Creating success for your patients and your practice.

The orthodontic business and treatment models are changing. Patients have more choice than ever, and they're playing an active role in their treatment decision making.

In the middle of all this change, our goals remain the same—creating beautiful, confident smiles that will last your patients a lifetime. We're here to support and empower your practice with everything you need:

- Educational materials to keep your practice ahead of the game.
- Hands-on training for you and your staff.
- Customizable marketing materials to drive more patient treatment.

Discover the difference at [3M.com/CeramicBrackets](https://www.3M.com/CeramicBrackets)



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-423-4588
[3M.com/Clarity](https://www.3M.com/Clarity)

3M, Clarity, and APC Adhesive
are trademarks of 3M.
Used under license in Canada.
© 3M 2020. All rights reserved.
70-2013-7245-0