

Technical information sheet

Solventum™ Microfluidic Diagnostic Tape 9793R

Effective: January 2010

Supersedes all previous information sheets

Product description

Microfluidic Diagnostic Tape 9793R consists of a clear polyolefin film coated on one side with a pressure sensitive acrylate adhesive. The tape is supplied on a white silicone coated paper release liner.

Applications

- Bioanalytical devices
- Diagnostic point of care test strips

Features and benefits

- High optical clarity
- Low auto-fluorescence
- “Hydrophobic” adhesive resists clouding
- Low level of extractables in a variety of solvents
- Adheres well to a wide range of substrates
- Suitable for use -70°C to 100°C (-94°F - 212°F)

Composition

Backing	2.0 mil (0.05 mm) clear polypropylene film
Adhesive	Pressure sensitive “Neutral” acrylic-based adhesive
Release liner	4.9 mil (0.12 mm) 63 lbs white silicone-coated paper

Typical product properties¹

Adhesion to steel²	25 ounces/inch width (7 N/25 mm width)
Liner release²	Maximum 50 g/in (0.5 N/25 mm)
Light transmission (Spectrophotometer)	92.3% at 700 nm (typical) 90.7% at 400 nm (typical)
Solvent extractables (Reference: GC/MS)	DMSO: None detected above background (typical) 80% ethanol/20% water: 5 µg/ml hydrocarbon acrylate esters (typical) 630 µg/ml antioxidant (typical)

¹ Test methods are available upon request.

² Lot specific results may be included on a Certificate of Analysis upon request at time of order.

Product dimensions: Product supplied on 3-inch (76 mm) white plastic cores & rolled liner side “in”.

Roll length	Maximum 500 yd (457 m) ± 2%
Roll width	Maximum 48 in (122 cm) ± 1/16 in (1.6 mm)

Note: For “EL” orders, product will be wound with the liner on the outside of the roll.

Splice Frequency: (Intended for continuous processing)

Roll length:	Maximum # of splices:
1 to 100 yd (0.9 to 91 m)	2
101 to 200 yd (92 to 183 m)	3
201 to 400 yd (184 to 366 m)	5
401 to 500 yd (367 to 457 m)	7

Please contact your Solventum Medical Surgical Account Representative for information on splice-free rolls.

Recommended storage & preliminary rollstock shelf life

Product as supplied by Solventum in original packaging will maintain stated test properties for a period of two (2) years from date stamped on shipping container when stored at temperatures of 50-80°F (10-27°C) and a relative humidity of 40-60 percent.

Please see reverse side for important product, safety and warranty information.

Product and safety information: User is solely responsible for determining the suitability of Solventum samples and products for the intended use including any necessary safety or toxicity assessment. Solventum will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

Warranty information

All statements, technical information and recommendations herein are based on tests Solventum believes to be reliable, but the accuracy or completeness thereof is not guaranteed. Solventum warrants only that products will meet Solventum's specifications at the time of shipment to the customer. Solventum does not offer any other warranty and does not warrant the performance, safety or such other characteristics of products in combination with other materials. Solventum specifically **does not** warrant products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of remedies: If products are proven not to meet Solventum's specifications, the sole and exclusive remedy available and Solventum's only obligation shall be, at Solventum's option, to replace such quantity of products which are proven out of specification or to refund the purchase price paid for products.

Limitations of liabilities: The remedies provided herein are exclusive remedies against Solventum for any alleged or actual nonconformance to specifications or defect or other failure in products or for Solventum's performance of its supply obligations. Under no circumstances is Solventum liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering information

Call us at Customer Service to place an order: 800-742-1994 (U.S.).

Visit our website: go.Solventum.com/MedTechOEM for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at Solventum Helpline 800-742-1994 (U.S.) or for international inquiries, please contact your local country representative.



Solventum MedTech OEM

Phone 1-800-742-1994

Web go.Solventum.com/MedTechOEM

©Solventum 2025. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates.

791981