



## 3M™ Sof-Lex™ Diamond Polishing System

# Polish with diamonds – skip the paste



How much time and effort do you spend creating beautiful smiles? Whether you currently use a rubberized finishing and polishing system or an intraoral diamond polish, the process can be time-consuming. And, even with your best effort, the gloss may not last. Solventum has a simple solution for both problems, using two of our innovative technologies.

### Restore with 3M™ Filtek™ Supreme Ultra Universal Restorative

Unsurpassed esthetics is just one reason why doctors use this nanocomposite. Thanks to Solventum's true nanotechnology, it is easy to polish and offers unsurpassed polish retention.

### Polish with the Sof-Lex Diamond Polishing System

Forget the messy paste. Our pre-polishing spiral prepares the restoration for final gloss, while our diamond-impregnated polishing spiral gives your restorations that gorgeous paste-like gloss. The system offers the convenience of a rubberized system while also adapting to all tooth surfaces.

You'll be happy to know that while the spirals are effective, they're also kinder to gingival tissues\*—and maintain the integrity and anatomy of your restorations!

When patients leave your office smiling, you'll marvel at how simple it's become to create beautiful, natural-looking esthetics.

\*Compared to other finishing and polishing tools.

## You can create a diamond paste-like gloss with just two steps



1

3M™ Sof-Lex™  
Pre-Polishing  
Spiral (beige)



2

3M™ Sof-Lex™  
Diamond Polishing  
Spiral (pink)

## A difference you can see!



**Filtek Supreme Ultra Universal Restorative** polished with **Sof-Lex Diamond Polishing System** (left) vs. TPH Spectra® Universal Composite polished with Enhance® Finishing System and PoGo® Polishing System (right). Notice a clearer reflection with the Sof-Lex Diamond Polishing System.

## Summary of advantages

- Imparts paste-like gloss in the convenience of a rubberized system
- Unique, flexible shape adapts to all tooth surfaces
- Fast and easy to use
- Multi-use, can be sterilized and reused
- High, long-lasting gloss when used with Filtek Supreme Ultra Universal Restorative

# Don't you love it when you find a product that works great, but is also cost effective?

Then you'll love this! The multipurpose 3M™ Sof-Lex™ spirals can be sterilized and reused, making them an economical choice for your practice.

Because abrasive particles are embedded throughout the rubberized spirals, all surfaces of the spiral design can be used: side, top, bottom and edge.

And forget the need to stop and switch between traditional points, cups, discs and brushes, too. The flexible spirals adapt to all tooth surfaces— anterior and posterior, convex and concave.



## 3M™ Sof-Lex™ Pre-Polishing Spiral

The pre-polishing spiral (beige) smooths and removes scratches in restorations that develop during contouring and prepares the surface for high-gloss polishing.



## 3M™ Sof-Lex™ Diamond Polishing Spiral

The diamond polishing spiral (pink) contains diamond particles to create a smooth, high-gloss polish.

**You're going to love how easy it is to create a beautiful smile!**



## Ordering information

Item number	Product information
5092-IM	<b>Sof-Lex Diamond Polishing System—Introductory kit</b> Kit includes: 5 Sof-Lex™ Pre-Polishing Spirals (Beige); 5 Sof-Lex™ Diamond Polishing Spirals (Pink); 1 Sof-Lex™ RA Mandrel; Instructions for Use; Technique card
5090	<b>Sof-Lex Pre-Polishing Spirals—Refill</b> Includes: 15 Sof-Lex Pre-Polishing Spirals (Beige)
5091	<b>Sof-Lex Diamond Polishing Spirals—Refill</b> Includes: 15 Sof-Lex Diamond Polishing Spirals (Pink)
1983RA	<b>3M™ Sof-Lex™ Mandrel—Refill (3)</b>

