



# Designed for mobile patients

## Solventum™ ActiV.A.C.™ Therapy System

Now including the first-ever Solventum™ V.A.C.®  
Therapy Dressing with up to 7-days of wear time!

# The Solventum™ ActiV.A.C.™ Therapy Unit is a portable negative pressure wound therapy device designed for the mobile patient

The features of this ergonomically designed therapy unit include:

## For the patient

- Lightweight. Weighs only ~1.21 kg
- Small size with a low profile that can be worn close to the body
- Easy-to-use, single-touch therapy on/off operation
- Alarm notifications that are easy to recognize and correct
- Easy, quick release 300 ml canister
- ~13 hour average battery for activities of daily living

## For the clinician

### Solventum™ SensaT.R.A.C.™ technology

Only Solventum™ V.A.C.® Therapy Systems provide patented SensaT.R.A.C. Technology, a real time pressure feedback system.

- Adjusts pump output to the wound site, compensating for wound distance, wound position, exudate characteristics, and patient movement
- Delivers and maintains programmed negative pressure to support optimal healing outcomes

### Solventum™ Seal Check™ Feature

Designed to help clinicians and patients identify and troubleshoot negative pressure leaks.

- Instant feedback. Audio and visual cues are delivered in real-time, allowing easy location of dressing leaks
- Enhanced therapy confidence. A proper seal helps maintain programmed therapy to support optimal wound healing
- Possible reduction of unscheduled visits. Easier for clinicians to help troubleshoot problems over the phone



## Settings guide

Recommended therapy settings and pressure ranges are pre-programmed for indicated wound types:

- Easy set up. With pre-programmed settings by wound type, clinicians can initiate therapy quickly and efficiently
- Ease of use. Intuitive control functions make it easier to use, especially for clinicians who are not as familiar with V.A.C.® Therapy

## Therapy history reporting

Document Patient Compliance with V.A.C.® Therapy.

- Monitors therapy usage, settings, durations, and dressing and canister changes to help verify patients are receiving the prescribed course of therapy
- USB port for convenient data downloads

## Simplified touch screen

Enhanced, simplified control functions:

- Easy navigation. With full-color interface with intuitive touch screen controls
- Simplified patient mode. Designed to be easy to use and understand

# Solventum™ V.A.C.® Dressing options with Solventum™ SensaT.R.A.C.™ Technology



## Solventum™ V.A.C.® Peel and Place Dressing Kit

All-in-one design makes application easy

- The perforated, non-adherent layer allows for extended wear time up to 7-days<sup>1</sup> and fewer dressing changes
- Patented, all-in-one dressing allows for easy, time-saving application with less training
- The simple, versatile design of V.A.C.® Peel and Place Dressing is compatible with multiple Solventum™ V.A.C.® Therapy Units



## Solventum™ V.A.C.® Granufoam™ Bridge Dressings

Ideal for use with wounds in pressure sensitive areas that requiring bridging

- Preassembled components specifically designed to simplify dressing application
- Integrated bridge allows for SensaT.R.A.C. Pad placement away from the wound site
- Wicking layers help ensure intact skin stays dry
- Allows V.A.C.® Therapy to be used in conjunction with off-loading therapies



## Solventum™ V.A.C.® Granufoam™ Silver Dressings

More than just a silver dressing

- Combines the benefits of V.A.C.® Therapy and silver in one dressing, with no need for adjunct silver layers
- Specifically designed for use with V.A.C.® Therapy Systems



## Solventum™ V.A.C.® Simplace™ Dressings

Designed to simplify the V.A.C.® Dressing placement process

- Helps to create an environment that promotes healing by facilitating granulation tissue formation
- Fewer steps, easier and faster application
- Spiral cut foam is simple to size; no scissors necessary
- Design allows for easier bridging



## Solventum™ V.A.C.® Granufoam™ Dressings

Advanced wound dressings to assist granulation tissue formation and enhance exudate removal

- Can be easily trimmed to fit the contours of deep or irregularly shaped wounds
- Can be customized for bridging techniques when treating multiple wounds



## Solventum™ V.A.C.® White Foam Dressings

Comfortable and versatile wound dressings

- May be used over split thickness skin grafts
- Helps bolster flaps and grafts
- Higher tensile strength (than V.A.C.® Granufoam Dressing) allows for easy removal from tunnels and undermining
- Recommended in situations where hypergranulation responses are likely



## Solventum™ SensaT.R.A.C.™ Pad

Only SensaT.R.A.C. Technology actively measures, monitors, and manages the accurate delivery of the physician's prescribed therapy to the wound site

- Designed with patient comfort in mind
- Flexible pad material for easy application over body contours
- Discreet, low profile design



## Solventum™ V.A.C.® Dressings and Accessories ordering information

Part number	Description	Case quantity
<b>Up to 7-day wear time dressing kits for Solventum™ V.A.C.® Therapy*</b>		
EZ5SML/5	Solventum™ V.A.C.® Peel and Place Dressing Kit - Small	5
EZ5MED/5	Solventum™ V.A.C.® Peel and Place Dressing Kit - Medium	5
EZ5LRG/5	Solventum™ V.A.C.® Peel and Place Dressing Kit - Large	5
<b>Up to 48-72 hours, or no less than 3 times per week dressing kits for V.A.C.® Therapy</b>		
M8275051/10	Solventum™ V.A.C.® Granufoam™ Dressing Kit - Small	10
M8275051/5	Solventum™ V.A.C.® Granufoam™ Dressing Kit - Small	5
M8275052/10	Solventum™ V.A.C.® Granufoam™ Dressing Kit - Medium	10
M8275052/5	Solventum™ V.A.C.® Granufoam™ Dressing Kit - Medium	5
M8275053/10	Solventum™ V.A.C.® Granufoam™ Dressing Kit - Large	10
M8275053/5	Solventum™ V.A.C.® Granufoam™ Dressing Kit - Large	5
M8275065/5	Solventum™ V.A.C.® Granufoam™ Dressing Kit - X-Large	5
M8275098/10	Solventum™ V.A.C.® Granufoam™ Silver Dressing/Dressing Kit - Small	10
M8275098/5	Solventum™ V.A.C.® Granufoam™ Silver Dressing/Dressing Kit - Small	5
M8275096/10	Solventum™ V.A.C.® Granufoam™ Silver Dressing/Dressing Kit - Medium	10
M8275096/5	Solventum™ V.A.C.® Granufoam™ Silver Dressing/Dressing Kit - Medium	5
M8275099/10	Solventum™ V.A.C.® Granufoam™ Silver Dressing/Dressing Kit - Large	10
M8275099/5	Solventum™ V.A.C.® Granufoam™ Silver Dressing/Dressing Kit - Large	5
M8275045/5	Solventum™ V.A.C.® Simplace™ Ex Dressing Kit - Medium	5
M8275046/5	Solventum™ V.A.C.® Simplace™ Ex Dressing Kit - Small	5
M8275042/10	Solventum™ V.A.C.® Granufoam™ Bridge Dressing Kit	10
M8275042/5	Solventum™ V.A.C.® Granufoam™ Bridge Dressing Kit	5
M6275033/10	Solventum™ V.A.C.® White Foam Dressing, Foam Only - Small**	10
M6275034/10	Solventum™ V.A.C.® White Foam Dressing, Foam Only - Large**	10
M8275067/10	Solventum™ V.A.C.® White Foam Dressing Kit - Large**	10
M8275067/5	Solventum™ V.A.C.® White Foam Dressing Kit - Large**	5
M8275068/10	Solventum™ V.A.C.® White Foam Dressing Kit - Small**	10
M8275068/5	Solventum™ V.A.C.® White Foam Dressing Kit - Small**	5
<b>Accessories for V.A.C.® Therapy</b>		
M6275009/10	Solventum™ V.A.C.® Drape	10
M6275026/10	Solventum™ V.A.C.® Gel	10
M6275066/10	Solventum™ V.A.C.® Y-Connector	10
M6275066/5	Solventum™ V.A.C.® Y-Connector	5
M6275069/10	Solventum™ V.A.C.® Tubing Cap	10
M6275069/5	Solventum™ V.A.C.® Tubing Cap	5
<b>Solventum™ ActiV.A.C.™ Canisters</b>		
M8275058/10	Solventum™ ActiV.A.C.™ 300 ml Canister with Gel	10
M8275058/5	Solventum™ ActiV.A.C.™ 300 ml Canister with Gel	5

\*See V.A.C.® Peel and Place Dressing IFU for more information on wounds appropriate for these dressings.

\*\* Less-adherent than V.A.C.® Granufoam Dressing



V.A.C.® Therapy has been awarded the APMA Seal of Acceptance, which recognizes products that have been found to promote quality foot health and be of significant value when used in regular professional treatment.



### Solventum Advanced Wound Care

12930 IH10 W  
San Antonio, TX 78249  
USA

Phone 800-275-4524

Web Solventum.com

#### Reference

1. SAT-BSER-05-869347 VAC Peel and Place BSER

**Note:** Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals. Rx Only.

**Note:** Disposable components of the Solventum™ V.A.C.® (Vacuum Assisted Closure®) Therapy System, including the foam dressing (i.e., Solventum™ V.A.C.® Granufoam™ Dressing, V.A.C.® Granufoam Silver Dressing, or V.A.C.® White Foam Dressing), tubing and drape are packaged sterile and are manufactured without natural rubber latex. Solventum™ V.A.C.® Therapy Unit canisters are packaged fluid path sterile and are manufactured without natural rubber latex. All disposable components of the V.A.C.® Therapy System are for single use only. To help ensure safe and effective use, the V.A.C.® Granufoam, V.A.C.® Granufoam Silver and V.A.C.® White Foam Dressings are to be used only with V.A.C.® Therapy Units.

© Solventum 2025. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates. Other trademarks are the property of their respective owners.

US\_70-2011-8287-3