



Designed to impress. Engineered for the practice.


Innovative orthodontic solutions to meet your needs.

CLARITY™ ALIGNERS
CLARITY™ PRECISION GRIP ATTACHMENTS
CLARITY™ DIGITAL BONDING
CLARITY™ RETAINERS
CLARITY™ PORTAL



Clarity™ Aligners

Customized clear aligner treatment to meet your needs and the needs of your patients




	Clarity Aligners Flex	Clarity Aligners Force
Film thickness	0.625 mm (thinner yet durable!)	0.75 mm
Max velocity	0.25 mm per stage 0.1 mm per overcorrection tray for space closure	
Force persistence	Excellent force persistence, maintaining a consistent force over two weeks.	

TX Design	Start with Clarity Aligners Flex	Start with Clarity Aligners Force
Tx Philosophy	Wire sequencing	Segmental mechanics
Patient malocclusion	Crowding, Undercuts	Deep curve of spee
Patient behavior	Clenching, Grinding	
Transition proposed by Solventum based on near completion of key movements	Rotation, Proclination, Extrusion	Expansion, Torque, Sequential bodily movement, Segmental intrusion


Clarity™ Precision Grip Attachments

Say goodbye to imprecision and hello to easier attachments

Clarity Precision Grip Attachments



Conventional Attachments



No composite flash removal



Reduce doctor time at the bonding appointment **by 49%**



Experience a reduction in treatment time



Attachments are pre-loaded – **saving 12 minutes per bonding procedure on average**




No composite flash removal – reducing unwanted tooth movements



Excellent stain resistance for clear, esthetic smiles during treatment




Impressive durability and high fracture strength helps the attachments maintain shape integrity



See the difference Clarity Precision Grip Attachments make:

The majority of evaluating doctors agree - Clarity Precision Grip Attachments can shorten treatment time for similar malocclusions vs. conventional attachments.



Clarity Precision Grip Attachments offer predictable engagement with aligners due to no composite flash.

Conventional attachments

- ✗ Unpredictable shapes
- ✗ Time intensive
- ✗ Unreliable aligner engagement
- ✗ Excess composite (flash)

Clarity Precision Grip Attachments

- ✓ Precisely customizable
- ✓ Simplified process
- ✓ Predictable aligner engagement
- ✓ No excess composite (flash)

Clarity™ Digital Bonding

Bond an entire arch of brackets simultaneously with this digitally designed custom device

APC™ Flash-Free Adhesive Coated Appliances offer resin-soaked nonwoven mats that adhere to enamel for a more reliable bond that also eliminates the process of flash removal.



The tray

Compatible with:



Clarity™ Advanced Ceramic Brackets



Victory Series™ Low Profile Brackets




SmartClip™ SL3 Self-Ligating Brackets

Made of two layers:

- Soft gel inlay to conform to the patient's anatomy, making removal easier
- Hard outer shell to add rigidity and support placement of the tray

Full arch tray allows for simplified bonding procedures where all brackets are placed at once



Blue midline marker to assist with tray alignment when seating

Soft inner tray does not bond with cured adhesive, minimizing the risk of bond failure

Clarity™ Digital Bonding supports bite blocks and bite ramps, giving you more flexibility in treatment planning. These features empower you to bond a full arch of brackets and bite features in a single procedure.

Doctors can expect to save treatment time with Clarity Digital Bonding (up to 6 months).

Clarity™ Retainers

The key to keeping a perfect smile?
Retention, clearly.

- Digitally remove appliances to retain fixed appliance and aligner cases
- Designed for optimal patient comfort
- Direct-to-patient delivery
- An all-digital approach – retainers can be ordered from a new scan or final aligner stage
- Increased durability to withstand more insertions and removals

All retainer orders come with a one-year warranty

Time Savings:

Doctor approval is not required for retainer orders, which accelerates the ordering and shipping process and reduces a step for your team.

“I wear Clarity Retainers, and the plastic is extremely comfortable. My patients are also thrilled with the lightweight feel and have never been more pleased.”

- Dr. David Pechersky
Pechersky Orthodontics

Clarity™ Portal

One integrated workflow



Cloud-based software provides convenient access to Tx Design



Accepts any high quality STL file



Integrated workflows supporting Tx Designs for Clarity™ Aligners with Clarity™ Precision Grip Attachments, Clarity™ Retainers and Clarity™ Digital Bonding



We're here for you

We're committed to your clinical and professional success. Reach out to your sales representative, clinical specialist, or call Clinician Support at **1-800-423-4588**.



Solventum Dental Solutions
1-800-423-4588
Solventum.com

© Solventum 2026. Solventum, the S logo, APC and Clarity are trademarks of Solventum or its affiliates. Other trademarks are the property of their respective owners. 70-2013-7338-1