



Littmann

Stethoscopes

Listen with precision

Littmann® Cardiology IV™ Stethoscope

Designed to help you listen for hard-to-hear sounds in critical environments, including those related to aortic stenosis, S3 gallop murmurs and faint pulmonary anomalies.



Available in 30 colours



Isolation of subtle sounds

Isolate subtle sounds critical for patient assessment in challenging environments, including those related to aortic stenosis, S3 gallop murmurs and faint pulmonary anomalies.



Versatile application

Hear sounds related to patient status changes in critical and dynamic care settings like the ED, ICU, Cardiac ICU and step-down units, empowering confident diagnoses and informed decision-making.



Enhanced design and comfort

Experience an ergonomic design with 40% larger chestpiece and 60% deeper bell than the Littmann® Classic III™ Stethoscope, suitable for both adult and paediatric patients.

Details

Acoustic performance	9 (1–10) scale
Chestpiece design/technology	Machined stainless steel, two-sided
Headset material	Lightweight aluminium alloy
Length/total weight	69cm/177g
Diaphragm details/diameter	Two, single-piece turnable diaphragms 4.3cm, 3.4cm
Warranty	7 years

The difference is in the details

Hear the hard stuff

Littmann® Cardiology IV™ Stethoscopes are more than twice as loud as the next leading stethoscope, allowing you to more reliably and consistently hear heart sounds at lower frequencies (below 120Hz) including those related to Korotkoff sounds, Mitral Stenosis and S3 and S4 Gallops when compared to other leading stethoscopes.*

Adjustable headset

Headset tension can be easily adjusted for individual fit and comfort. Angled eartubes align with ear canals for comfort and maximum sound transmission.

Lightweight design

Littmann stethoscopes are lighter than other stethoscopes** and designed to feel comfortable when draped around your neck or carried in a pocket, giving you an extra step in your day.

Tunable diaphragm technology

Hear different frequency sounds without turning over and repositioning the chestpiece. Just hold the chestpiece with light pressure to hear low frequency sounds, or press a little more to hear higher frequency sounds.

Handcrafted quality

Littmann® Stethoscopes are designed and manufactured in the United States with globally sourced materials following rigorous processes and strict controls. The result? Littmann stethoscopes are reliably louder than other stethoscopes.*

Precision chestpiece

Designed to provide excellent sound quality over a lifetime of use, chestpieces are crafted with stainless steel and high-density zinc alloy due to their excellent acoustic properties and their strength to resist impacts, scratches and chemicals. Double-sided chestpieces have precision stems that rotate smoothly and click into position.



Find yours at [Littmann.com/en-gb](https://www.littmann.com/en-gb)



**Authorised
Dealer**

*Based on tests against globally-available comparable stethoscopes in an equivalent class using recorded heart sounds with diaphragm.

**Based on published weights of globally-available comparable stethoscopes in an equivalent class.

© Solventum 2026. Solventum, the S logo, Littmann, the L Littmann logo, the L, Classic III and Cardiology IV are marks and/or registered marks of Solventum or its affiliates. Other trademarks are the property of their respective owners. OMG2302398.