

Littmann®

Stethoscopes

Choose the right stethoscope to help you perform at your best

	Littmann® CORE Digital Stethoscope	Littmann® Master Cardiology™ Stethoscope	Littmann® Cardiology IV™ Stethoscope	Littmann® Classic III™ Stethoscope	Littmann® Classic II™ Paediatric and Infant Stethoscopes	Littmann® Lightweight II S.E. Stethoscope
	<p>Switch from analogue listening to digital. Connects to Eko software on a smart device** to capture, store and share data.</p> <p>Chestpiece diameter double-sided Adult: 4.3cm Paediatric: 3.4cm</p> 	<p>Best acoustic response in Littmann mechanical line.</p> <p>Chestpiece diameter single-sided 4.8cm</p> 	<p>Diagnose and make decisions with confidence.</p> <p>Chestpiece diameter double-sided Adult: 4.3cm Paediatric: 3.4cm</p> 	<p>Monitor and assess a wide range of patients.</p> <p>Chestpiece diameter double-sided Adult: 4.3cm Paediatric: 3.4cm</p> 	<p>Specially designed for the smallest patients.</p> <p>Chestpiece diameter double-sided Paediatric: 3.4cm, 2.8cm Infant: 2.8cm, 1.7cm</p> 	<p>Teardrop-shaped chestpiece fits under blood pressure cuff.</p> <p>Chestpiece diameter double-sided 5.2cm x 4.5cm 3.5cm x 3.2cm</p> 
Acoustic performance level	10+	10	9	8	7	6
Approx. length	69cm	69cm	69cm	69cm	71cm	71cm
Approx. weight	232g	185g	177g	150g	105g/95g	118g
Warranty	2 years	7 years	7 years	5 years	3 years	2 years
Tunable diaphragm	●	●	●	●		●
Up to 40x amplification*	●					
Noise cancellation	●					
Mobile connectivity	●					
Dual-lumen tubing	●	●	●			
Additional special finishes	●	●	●	●	●	
Convertible open bell	●		●	●		
Soft sealing eartips	●	●	●	●	●	●
Dual-sided chestpiece	●		●	●	●	●

*At peak frequency vs. analogue mode.
**Subscription may be required.

© Solventum 2025. Solventum, the S logo and Littmann, the L Littmann logo, the L, Master Cardiology, Cardiology IV and Classic III are trademarks of Solventum or its affiliates. Other trademarks are the property of their respective owners. OMG2019484.