

A full line of restorative solutions to fit your needs

3M™ Filtek™ Dental Restoratives have been used in more than 2 billion restorations worldwide. Our family of Filtek dental restoratives lets you choose just the right option to fit any clinical need – from highly aesthetic, to highly simple for an easier day.

	Aesthetics and versatility		Aesthetics and simplicity		Efficiency and simplicity	
Properties	3M™ Filtek™ Supreme XTE Universal Restorative	3M™ Filtek™ Supreme Flowable Restorative	3M™ Filtek™ Universal Restorative	3M™ Filtek™ Easy Match Universal Restorative	3M™ Filtek™ One Bulk Fill Restorative	3M™ Filtek™ Bulk Fill Flowable Restorative
 Primary benefits	• Beautiful aesthetics with natural-looking fluorescence • Outstanding strength for anterior and posterior	• Virtually no bubbles or 'run-on' • Ergonomic design • Selection of cannula tip sizes (20G, 27G)	• Universal opacity for consistent aesthetics • Pink Opaquer to mask dark areas • Low shrinkage stress	• Quick, easy and intuitive shade selection without sacrificing opacity • Simplified single-layer restorations, with no blocker needed • Naturally-adaptive opacity	• Increased opacity for improved aesthetics • One-step placement up to 5 mm	• Virtually no bubbles or 'run-on' • One-step placement up to 4 mm • Ergonomic design
 Number of shades	36 (Dentine, Body, Enamel, Translucent)	12	8 (plus Extra White and Pink Opaquer)	3	5	4
 TRUE nanotechnology	Yes	Yes	Yes	Yes	Yes	No
 Handling	Firm	Flows on demand	Creamy	Firm	Sculptable	Adapts easily
 Polish retention	Brilliant					N/A*
 Wear resistance	Excellent					N/A*
 Safe for composite warming and injection molding ^{1,2}	Capsules 	Syringes 	Capsules 	Capsules 	Capsules 	Syringes 

1. Warm capsules to up to 70°C for up to 1 hour. 2. Warm flowable syringes to up to 70°C for up to 1 hour, up to 25 times. 3. Internal data.

*Filtek Bulk Fill Flowable's intended use requires allowing at least 2 mm on the occlusal cavosurface for a universal or posterior composite.

We have you covered for every class of direct restoration

	Aesthetics and versatility		Aesthetics and simplicity		Efficiency and simplicity	
	3M™ Filtek™ Supreme XTE Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Supreme Flowable Restorative (Incremental up to 2 mm)	3M™ Filtek™ Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Easy Match Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)	3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I 	✓	✓	✓	✓	✓ (up to 4 mm)	✓ (Liner/Base)
Class II 	✓	✓	✓	✓	✓ (up to 5 mm)	✓ (Liner/Base)
Class III 	✓	✓	✓	✓		✓
Class IV 	✓	✓	✓	✓		
Class V 	✓	✓	✓	✓		✓
Minimally invasive	✓	✓	✓	✓	✓	✓
Stress bearing	✓	✓	✓	✓	✓	

✓ = recommended use



Solventum Dental Solutions

KCI Medical Australia Pty Ltd

Level 3, Building A, 1 Rivett North Ryde, NSW 2113 | 1300 363 454

KCI New Zealand Unlimited

Suite 1701, Level 17, PwC Tower, 15 Customs Street West
Auckland Central, Auckland 1010 | 0800 808 182

www.solventum.com

Refer to Instructions for Use for full list of indications.



Request a demo

© Solventum 2026. Solventum, the S logo and Filtek are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. All other trademarks are owned by their respective owners.