

A full line of restorative solutions to fit your needs

3M™ Filtek™ Dental Restoratives have been used in more than 2 billion restorations worldwide. Our family of Filtek dental restoratives lets you choose just the right option to fit any clinical need – from highly aesthetic, to highly simple for an easier day.

		Aesthetics and versatility		Aesthetics and simplicity		Efficiency and simplicity	
Properties		3M™ Filtek™ Supreme XTE Universal Restorative	3M™ Filtek™ Supreme Flowable Restorative	3M [™] Filtek [™] Universal Restorative	3M™ Filtek™ Easy Match Universal Restorative	3M™ Filtek™ One Bulk Fill Restorative	3M™ Filtek™ Bulk Fill Flowable Restorative
	Primary benefits	Beautiful aesthetics with natural-looking fluorescence Outstanding strength for anterior and posterior	 Virtually no bubbles or 'run-on' Ergonomic design Selection of cannula tip sizes (20G, 27G) 	Universal opacity for consistent aesthetics Pink Opaquer to mask dark areas Low shrinkage stress	Quick, easy and intuitive shade selection without sacrificing opacity Simplified single-layer restorations, with no blocker needed Naturally-adaptive opacity	Increased opacity for improved aesthetics One-step placement up to 5 mm	Virtually no bubbles or 'run-on' One-step placement up to 4 mm Ergonomic design
F	Number of shades	36 (Dentine, Body, Enamel, Translucent)	12	8 (plus Extra White and Pink Opaquer)	3	5	4
	TRUE nanotechnology	Yes	Yes	Yes	Yes	Yes	No
	Handling	Firm	Flows on demand	Creamy	Firm	Sculptable	Adapts easily
₩	Polish retention			Brilliant			N/A*
(S)	Wear resistance			Excellent			N/A*
I E	Safe for composite warming and injection molding ^{1,2}	Capsules	Syringes	Capsules	Capsules	Capsules	Syringes

^{1.} Warm capsules to up to 70°C for up to 1 hour. 2. Warm flowable syringes to up to 70°C for up to 1 hour, up to 25 times. 3. Internal data.

^{*}Filtek Bulk Fill Flowable's intended use requires allowing at least 2 mm on the occlusal cavosurface for a universal or posterior composite.

We have you covered for every class of direct restoration

	Aesthetics and versatility		Aesthetics and simplicity		Efficiency and simplicity	
	3M [™] Filtek [™] Supreme XTE Universal Restorative (Incremental up to 2 mm)	3M [™] Filtek [™] Supreme Flowable Restorative (Incremental up to 2 mm)	3M [™] Filtek [™] Universal Restorative (Incremental up to 2 mm)	3M [™] Filtek [™] Easy Match Universal Restorative (Incremental up to 2 mm)	3M [™] Filtek [™] One Bulk Fill Restorative (Bulk up to 5 mm)	3M [™] Filtek [™] Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I	\bigcirc	\bigcirc	\bigcirc	\bigcirc	(up to 4 mm)	(Liner/Base)
Class II	\bigcirc	\bigcirc	\bigcirc	\bigcirc	(up to 5 mm)	(Liner/Base)
Class III	\bigcirc	\bigcirc	\bigcirc	\bigcirc		\bigcirc
Class Ⅳ	\bigcirc	\bigcirc	\bigcirc	\bigcirc		
Class ☑	\bigcirc	\bigcirc	\bigcirc	\bigcirc		\bigcirc
Minimally invasive	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc
Stress bearing	\bigcirc	\bigcirc	\bigcirc	\bigcirc	\bigcirc	






Solventum Dental Solutions

KCI Medical Australia Pty Ltd Level 3, Building A, 1 Rivett North Ryde, NSW 2113 | 1300 363 454

KCI New Zealand Unlimited Suite 1701, Level 17, PwC Tower, 15 Customs Street West Auckland Central, Auckland 1010 | 0800 808 182 Refer to Instructions for Use for full list of indications.





© Solventum 2025. Solventum, the S logo and Filtek are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. All other trademarks are owned by their respective owners.