






























## ECG electrodes selection guide

# Using the right electrode matters

Better patient monitoring begins by choosing the right ECG electrode. When you choose from the Solventum Red Dot ECG Electrodes portfolio, you'll find the right electrode for a variety of patients, skin types, and monitoring situations. As a leader in adhesive innovation, Solventum knows how to make Red Dot electrodes stick and stay in place. The quality of patient monitoring directly affects the quality of patient care. Red Dot electrodes are easy to apply and provide uninterrupted patient monitoring with a trace quality you can trust.



	Adult						Adult/Pediatric	Pediatric	Neonatal	
	Universal							Diagnostic		
	2231	2228	2560 (2560, 2560-3, 2560-5)	2660 (2660-3, 2660-5)	2238	2289 2289T, 2289PA, 2289PAL	2268 (2268-3, 2268-5)	2330	2248	2269T
										
										
	Diaphoretic	Diaphoretic	Diaphoretic	Repositionable Fragile skin (2660 only)*	Fragile skin*		Fragile skin*	Fragile skin	Pediatric	N/A
<b>Adhesion</b>		 Same as 2228/2244		<b>2660:</b>  <b>2670:</b> 			 Gentle	<b>2330:</b>  <b>2360:</b> 	 Same adhesive as 2239 but smaller footprint	
<b>Trace quality</b>	<b>2231:</b> ●●●●●● <b>2244:</b> ●●●●●●	●●●●●●	●●●●●●	●●●●●●	●●●●●●	●●●●●● Same adhesive as 2670	●●●●●●	●●●●●●	Use a gentler solid gel	●●●●●●
<b>MR conditional</b>	<b>2231:</b> No <b>2244:</b> Yes	No	<b>2560:</b> No <b>2570:</b> Yes	Yes	No	Yes - Not sure if tested?	Yes	Yes	No	Yes
<b>Radiolucent</b>	<b>2231:</b> No <b>2244:</b> Yes	No	<b>2560:</b> No <b>2570:</b> Yes	Yes	No	Yes - Not sure if tested?	Yes	Yes	No	Yes
<b>Backing</b>	<b>2231:</b> Soft cloth <b>2244:</b> Foam	Foam	Foam	Soft cloth	<b>2238:</b> Soft cloth <b>2239:</b> 3M™ Micropore™ Tape	Clear tape	Soft cloth	Polyester film	Soft cloth	Clear tape
<b>Bag quantity</b>	50	<b>2228:</b> 50 <b>2228B:</b> 40	3, 5, 50	3, 5	<b>2238:</b> 50 <b>2239:</b> 50	<b>2289T:</b> 3 <b>2289PA:</b> 5 <b>2289PAL:</b> 5	3, 5	100	50	3
<b>Opened pouch expiration</b>	30 days	30 days	45 days	Use immediately	<b>2238:</b> 30 days <b>2239:</b> 30 days	Use immediately	Use immediately	7 days	30 days	Use immediately

© Solventum 2026. Solventum and the S logo are trademarks of Solventum or its affiliates. 3M and the 3M logo are licensed trademarks of 3M.