

 solventum

3M™ RelyX™ Universal
Resin Cement

Goodbye
Complexity



The two that can do it all

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive.

Together they eliminate the hassle and confusion of multiple resin cements, primers and adhesives, simplifying your direct and indirect restorative workflows. With fewer products to stock and replenish, you save space, time and money.

So simple yet so powerful

- Two components for virtually all indirect and direct bonding indications
- Outstanding self-adhesion of RelyX Universal Resin Cement
- Scotchbond Universal Plus Adhesive further enhances bond strength to all substrates
- Unique 3M automix syringe
- Easy excess clean-up
- High esthetics



Less stock, more savings

95% of dentists agreed the universal character of RelyX Universal Resin Cement helps to reduce stock and save costs.

Game-changing syringe

Proprietary 3M design takes automix cement applications to the next level. Developed with your needs in mind for easy dispensing, cleaning and storage. The small, low-waste 3.4g syringe provides enough material for around 15 applications.

Easy excess clean-up

Efficient to the last step. It's unique chemistry ensures excess outflow stays put at the restoration margins, and can be easily removed after tack curing.



Self-adhesive cementation of a 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia bridge; Image courtesy of Dr. Jan Gueth, Germany



Easier for operator staff

97% of dentists agreed using one universal versus multiple resin cements is easier for operator staff.

“Excellent hygienic dispensing system with little waste.”
M.W., Great Britain***

“Easy handling, few components and therefore less potential for error.”
J.L., Germany ***

80%

less cement waste* with Micro Mixing Tip

More hygienic.**
Self-sealing for storage without used mixing tip.

More ergonomic shape with 50% less plastic waste*

Thin, long elongation tip for easy root canal application

The strong partner for adhesive cases



3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque, all-in-one universal adhesive for all direct and indirect bonding indications. And, it is a fully aligned system: the adhesive is cured by the cement, so no light-cure is needed.



Scotchbond Universal Plus Adhesive serves as universal primer for all materials including glass ceramics and further boosts the cement's excellent self-adhesive bond e.g. to zirconia.

Outstanding self-adhesion for maximized efficiency


With its superior self-adhesive bond strength to dentin, RelyX Universal Resin Cement reliably treats most cases with no need for an additional adhesive or primer. For light- or self-cure.



Virtually no post-operative sensitivity

99% of dentists reported that RelyX Universal Resin Cement shows virtually no post-operative sensitivity***

High esthetics with four color-stable, fluorescent shades



Ordering information

3M™ RelyX™ Universal Resin Cement

Item no. Product information

Intro Kit

56968 4x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60x 3M™ RelyX™ Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

Trial Kits

Each trial kit contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

56969 **Shade TR** **56970** **Shade A1**

Refills

Each refill contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips

56971 **Shade TR** **56973** **Shade A3O**

56972 **Shade A1** **56974** **Shade WO**

Value Packs

Each value pack contains: 3x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45x 3M™ RelyX™ Universal Micro Mixing Tips, 15 Elongation Tips

56977 **Shade TR** **56978** **Shade A1**

Mixing Tips

56975 30x 3M™ RelyX™ Universal Micro Mixing Tips

56976 15x 3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips

3M™ Scotchbond™ Universal Plus Adhesive

Vial

41293 **Intro Kit:** 1x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips

41294 **Refill:** 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)
41294 **1x 5 ml vial** **41295** **3x 5 ml vial**

Unit Dose

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose, 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips

41297 **Intro Kit: 50x L-Pop™** **41299** **200x L-Pop™**

41298 **100x L-Pop™** **41304** **400x L-Pop™**

3M™ Scotchbond™ Universal Etchant Etching Gel

41263 2x 3M™ Scotchbond™ Universal Etchant (3 ml each), 50 dispensing tips



3M™ RelyX™ Universal Resin Cement Intro Kit (56968)



3M™ RelyX™ Universal Resin Cement Trial Kit TR (56969)



3M™ RelyX™ Universal Resin Cement Refill TR (56971)



3M™ Scotchbond™ Universal Plus Adhesive Refill (41294)



3M™ Scotchbond™ Universal Plus Adhesive Unit Dose:

For efficient hygiene management. Reduces the risk of cross-contamination.

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose 100 L-Pop (41298)



Solventum Dental Solutions
 1-800-634-2249
Solventum.com

For more information, visit Solventum.com

© Solventum 2026. Solventum, the S logo, L-Pop and RelyX are trademarks of Solventum or its affiliates. 3M, the 3M logo and Scotchbond are trademarks of 3M. All other trademarks are owned by their respective owners.