


# 3M™ Tube Securement Device

Strong hold. Gentle Care.

Designed to help prevent pressure injuries and medical adhesive-related skin injury (MARSI), this strong, easy-to-use solution secures and stabilizes non-critical IV tubes and fluid drainage tubes made from silicone, PVC and silicone coated latex spanning 6-24 French.



**Skin performance:** *The feather designates products that deliver the securement power you need while minimizing damage to skin.*

## Reliable securement

Formulated with a strong but gentle adhesive to secure medical tubes and catheters and be changed as clinically indicated.

## Easy to use

Easy to apply and removes cleanly in one piece; allows for standardization of technique.

## Skin safe

Non-sensitizing silicone adhesive has low potential to cause dermatitis with minimal skin damage and pain upon removal.<sup>1</sup>

## Helps reduce risk

With no hard components to cause pressure, product's shape minimizes contact between skin and tube reducing risk of pressure injury.

# Indications and Applications

When it comes to securing lines tubes and drains, there are a variety of considerations, including securement performance, ease of use, standardization of practice, and skin safety.



## Vascular access

- The tube securement device meets the Infusion Therapy Standards of Practice guidelines as an adhesive securement device and should be used in conjunction with sutures and/or primary securement dressings.
- Consider using the securement device as secondary securement for primary tubing and extension tubing on central lines and peripheral lines.
- It can help organize lines with catheters that have multiple lumens.

## Surgical drains

- Provides additional securement in conjunction with sutures to fluid drainage tubing to help prevent movement.
- Can provide secondary securement for negative pressure wound therapy drain tubes.



## Urinary catheters

- 9x faster removal.<sup>2</sup>
- No removal agent required.

To order, contact your Solventum Sales Representative or call 1-800-228-3957.

Refer to product Instructions for Use (IFU) for other important information.

## Ordering information

Product	Product name	Product number	Tube and catheter size	Product size	Devices/box	Boxes/case
	3M™ Tube Securement Device, small	1500U	6-12 French	1.7 in x 2.4 in (4.3 cm x 6.1 cm)	25	4
	3M™ Tube Securement Device, large	1501U	12-24 French	2.2 in x 3.5 in (5.5 cm x 9.0 cm)	25	4

2. When compared to StatLock® Foley Stabilization Device in a comparative prospective study of 24 healthy human subjects.

